
# Aesculap<sup>®</sup> Vitelene<sup>®</sup>

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- Kurtz SM, Gawel HA, Patel JD.  
History and Systematic Review of Wear and Osteolysis Outcomes for First-generation Highly Crosslinked Polyethylene.  
*Clin Orthop Relat Res.* 2011 Aug;469(8):2262-77.
- Jäger M, Behringer M, Zilkens C, Matheney T, Krauspe R.  
Initial increased wear debris of XLPE-Al2O3 bearing in total hip arthroplasties.  
*Arch Orthop Trauma Surg.* 2010 Dec;130(12):1481-6.
- Whittaker JP, Charron KD, McCalden RW, Macdonald SJ, Bourne RB.  
Comparison of steady state femoral head penetration rates between two highly cross-linked polyethylenes in total hip arthroplasty.  
*J Arthroplasty.* 2010;25:680-6.
- Calvert GT, Devane PA, Fielden J, Adams K, Horne JG.  
A double-blind, prospective, randomized controlled trial comparing highly cross-linked and conventional polyethylene in primary total hip arthroplasty.  
*J Arthroplasty.* 2009;24:505-10.
- Garvin KL, Hartman CW, Mangla J, Murdoch N, Martell JM.  
Wear analysis in THA utilizing oxidized zirconium and crosslinked polyethylene.  
*Clin Orthop Relat Res.* 2009;467:141-5.
- Garcia-Rey E, Garcia-Cimbrelo E, Cruz-Pardos A, Ortega-Chamarro J.  
New polyethylenes in total hip replacement: a prospective, comparative clinical study of two types of liner.  
*J Bone Joint Surg Br.* 2008;90:149-53.
- Glyn-Jones S, McLardy-Smith P, Gill HS, Murray DW.  
The creep and wear of highly cross-linked polyethylene: a three-year randomised, controlled trial using radiostereometric analysis.  
*J Bone Joint Surg Br.* 2008;90:556-61.
- Bitsch RG, Loidolt T, Heisel C, Ball S, Schmalzried TP.  
Reduction of osteolysis with use of Marathon cross-linked polyethylene: a concise follow-up, at a minimum of five years of a previous report.  
*J Bone Joint Surg Am.* 2008;90:1487-91.
- Bragdon CR, Kwon YM, Geller JA, Greene ME, Freiberg AA, Harris WH, Malchau H.  
Minimum 6-year followup of highly cross-linked polyethylene in THA.  
*Clin Orthop Relat Res.* 2007;465:122-7.
- Digas G, Karrholm J, Thanner J, Herberts P.  
5-year experience of highly cross-linked polyethylene in cemented and uncemented sockets: two randomized studies using radiostereometric analysis.  
*Acta Orthop.* 2007;78:746-54.

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- Oral E, Ghali BW, Neils A, Muratoglu OK.  
A new mechanism of oxidation in ultrahigh molecular weight polyethylene caused by squalene absorption.  
*J Biomed Mater Res B Appl Biomater.* 2012 Apr;100(3):742-51. doi: 10.1002/jbm.b.32507.
- Micheli BR, Wannomae KK, Lozynsky AJ, Christensen SD, Muratoglu OK.  
Knee simulator wear of vitamin e stabilized irradiated ultrahigh molecular weight polyethylene.  
*J Arthroplasty.* 2012 Jan;27(1):95-104.
- Bracco P, Oral E.  
Vitamin E-stabilized UHMWPE for total joint implants: a review.  
*Clin Orthop Relat Res.* 2011 Aug;469(8):2286-93.
- Rowell S, Oral E, Muratoglu O.  
Comparative Oxidative Stability of  $\alpha$ -Tocopherol Blended and Diffused UHMWPE at 3 Years of Real-Time Aging.  
*J Orthop Res.* 2011 May;29(5):773-80. doi: 10.1002/jor.21288.
- Oral E, Muratoglu OK.  
Vitamin E diffused, highly crosslinked UHMWPE: a review.  
*Int Orthop.* 2011 Feb;35(2):215-23.
- Currier BH, Van Citters DW, Currier JH, Collier JP.  
In vivo oxidation in remelted highly cross-linked retrievals.  
*J Bone Joint Surg Am.* 2010 Oct 20;92(14):2409-18.
- Gomez-Barrena E, Puertolas JA, Munuera L, Kontinen YT.  
Update on UHMWPE research: from the bench to the bedside.  
*Acta Orthop.* 2008 Dec;79(6):832-40.
- Dumbleton JH, D'Antonio JA, Manley MT, Capello WN, Wang A.  
The basis for a second-generation highly cross-linked UHMWPE.  
*Clin Orthop Relat Res.* 2006 Dec;453:265-71.
- Wannomae KK, Christensen SD, Freiberg AA, Bhattacharyya S, Harris WH, Muratoglu OK.  
The effect of real-time aging on the oxidation and wear of highly cross-linked UHMWPE acetabular liners.  
*Biomaterials.* 2006 Mar;27(9):1980-7.

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B. Braun Austria GmbH | Aesculap Division | Otto Braun-Straße 3-5 | 2344 Maria Enzersdorf  
Tel. +43 2236 4 65 41-0 | Fax +43 2236 4 65 41-177 | [www.bbraun.at](http://www.bbraun.at)

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B. Braun Medical AG | Aesculap Division | Seesatz 17 | 6204 Sempach  
Tel. +41 58 258 50 00 | Fax +41 58 258 60 00 | [www.bbraun.ch](http://www.bbraun.ch)

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Deutschland  
Tel. (0 74 61) 95-0 | Fax (0 74 61) 95-26 00 | [www.aesculap.de](http://www.aesculap.de)

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